

[Surgical Mask]

We provide array of medical masks. Anti-fog film, foam and tapes reduce distractions from fogging issues. Transparent and anti-scratching visor keeps eyes clear of blood and splash.

Material: PP+Meltblown+PP, 25+25+25gsm by request

Style: 3-Ply, ear-loop

Size: 17.5*9.5cm

Color:

Package: 50 pcs / box, 2000 pcs / ctn

Standard:

EN 14683 Type I, Type II, Type IIR

ASTM F2100 Level 1, Level 2, Level 3



REF	Features	Style	Color	Size
WLM2002	-	ear-loop	Blue	17.5*9.5cm
WLM2002S	Sponge	ear-loop	Blue	17.5*9.5cm