

### **[FFP3 Valved Respirator]**

#### **FFP3 NR D**

PFE:  $\geq 99\%$

Style: with valve

Standard: **EN149:2001+A1:2009**

Package: 10 pcs / box, 120 pcs / ctn

Style: metal nose clip, adjustable  
spandex ear-loop



**REF: WL-2012V-FFP3**

### **[FFP3 Valved Respirator with Acticarbon Layer]**

#### **FFP3 NR D**

PFE:  $\geq 99\%$

Style: with valve, with acticarbon layer

Standard: **EN149:2001+A1:2009**

Package: 10 pcs / box, 120 pcs / ctn

Style: metal nose clip, adjustable  
spandex ear-loop



**REF: WL-2012C-FFP3**