

WLG1005-L1 [Light Isolation Gown]

- Economic Impervious Fabric
- Fluid Repellent
- Low-linting
- Latex Free
- Acceptable Tearing Resistance
- Elastic Cuffs
- Non-sterile

Seams: Ultrasonic Welding

REF	Cuff	Color	Size
WLG1005-L1W	Elastic Cuff	White	S-3XL

Standard:

AAMI Level 1 Isolation Gown



WLG1005-L2 [Standard Isolation Gown]

- Impervious Fabric
- Fluid Repellent
- Low-linting
- Latex Free
- Ideal Tearing Resistance
- Knitted Cuffs / Elastic Cuffs
- Non-sterile

Seams: Ultrasonic Welding

REF	Cuff	Color	Size
WLG1005-L2KCY	Knitted Cuff	Yellow	S-3XL

Standard:

AAMI Level 2 Isolation Gown



WLG1005-L4 [Heavy Taped Isolation Gown]

- Heat Taping Technology
- Reinforced Seam Protection
- Heavy Duty Impervious Fabric
- Fluid Repellent
- Low-linting
- Latex Free
- Excellent Tearing Resistance
- Stitched and Taped
- Knitted Cuffs
- Non-sterile

Seams: Ultrasonic Welding

REF	Cuff	Color	Size
WLG1005-HiKCB	Knitted Cuff	Blue	S-3XL
WLG1005-L4KCB	Knitted Cuff	Blue	S-3XL

Standard:

AAMI Level 4 Isolation Gown
EN 13795 High Performance



WLG01 [PP Isolation Gown]

- Fabric: PP 20-60gsm
- Seam Type: Stitched
- Style: One piece opening at back, ties at waist
- With ties/velcro at back neck
- Knitted/elastic cuffs
- Color: ● ○ ● ●

Below are very basic referential styles, for additional requirements, please let us know.

REF	Weight	Style	Color	Size
WLG01B	20-60gsm	Elastic Cuff	Blue	S-3XL
WLG01KC	20-60gsm	Knitted Cuff	White	S-3XL



*All gowns are available with different colors and weights per your request.