

## [PP+PE Isolation Gown] WLG1005

Standard: AAMI PB70 Level 1, 2, 3, 4,

PP+PE fabric is born with the best water repellency and light weight. It is considered a no-brainer solution for protection against light liquid splash in most workplaces.

Fabric: PP+PE 30-60gsm

Seam Type: Stitched/ultrasonic welded

Style: One piece opening at back, ties at waist

With ties/velcro at back neck

Knitted/elastic cuffs

Color: ●●○○●

REF	Weight	Style	Color	Size
WLG1005	30-60gsm	Elastic Cuff	White	S-3XL
WLG1005B	30-60gsm	Elastic Cuff	Blue	S-3XL
WLG1005KC	30-60gsm	Knitted Cuff	White	S-3XL
WLG1005KCB	30-60gsm	Knitted Cuff	Blue	S-3XL
WLG1005KCY	30-60gsm	Knitted Cuff	Yellow	S-3XL
WLG1005Y	30-60gsm	Elastic Cuff	Yellow	S-3XL

