

## [Reinforced Convenience Gown] WLG06

One-step molding technology, heat sealing sewing without stitches

Over the head style to save time

The tear-away seam in the back for easy removal

Thumb hook and elastic cuff to prevent sleeve sliding

Chest and sleeves reinforced for an excellent barrier against fluids and pathogens.

Fabric: SMS/PP+PE 25gsm~45gsm

Size: 110x106cm, or customized

Color: 



REF	Fabric	Weight	Style	Color	Size
WLG06PECB	PP+PE Coated	25-55gsm	Elastic Cuff	Blue	110x106cm
WLG06PECY	PP+PE Coated	25-55gsm	Elastic Cuff	Yellow	110x106cm
WLG06PECTHB	PP+PE Coated	25-55gsm	Thumb Hoop	Blue	110x106cm
WLG06PECTHY	PP+PE Coated	25-55gsm	Thumb Hoop	Yellow	110x106cm
WLG06MB	SMS	25-55gsm	Elastic Cuff	Blue	110x106cm
WLG06MY	SMS	25-55gsm	Elastic Cuff	Yellow	110x106cm
WLG06MTHB	SMS	25-55gsm	Thumb Hoop	Blue	110x106cm
WLG06MTHY	SMS	25-55gsm	Thumb Hoop	Yellow	110x106cm