

WLG04

## 【Scrub Suit】

PP or SMS Fabric

Low-linting

Latex Free

Stitched

Elastic or tie waist pants

Short Sleeve With or Without Pocket

V Shape or Round Shape Collar

REF	Material	Weight	Style	Color	Size
WLG04MRB	SMS	30-40gsm	Round Collar	Blue	S-2XL
WLG04MVB	SMS	30-40gsm	V-shape Collar	Blue	S-2XL
WLG04PRB	PP	20-35gsm	Round Collar	Blue	S-2XL
WLG04PVB	PP	30-40gsm	V-shape Collar	Blue	S-2XL

\*All gowns are available with different colors and weights per your request.

