

[Standard Convenience Gown]

WLG05

One-step molding technology, heat sealing sewing without stitches

Over the head style to save time

The tear-away seam in the back for easy removal

Thumb hook and elastic cuff to prevent sleeve sliding

Fabric: SMS/PP+PE 25gsm~45gsm

Size: 110x106cm, or customized

Color: ●●



REF	Fabric	Weight	Style	Color	Size
WLG05MB	SMS	25-55gsm	Elastic Cuff	Blue	110x106cm
WLG05MTHB	SMS	25-55gsm	Thumb Hoop	Blue	110x106cm
WLG05MTHY	SMS	25-55gsm	Thumb Hoop	Yellow	110x106cm
WLG05MY	SMS	25-55gsm	Elastic Cuff	Yellow	110x106cm
WLG05PECB	PP+PE Coated	25-55gsm	Elastic Cuff	Blue	110x106cm
WLG05PECTHB	PP+PE Coated	25-55gsm	Thumb Hoop	Blue	110x106cm
WLG05PECTHY	PP+PE Coated	25-55gsm	Thumb Hoop	Yellow	110x106cm
WLG05PECY	PP+PE Coated	25-55gsm	Elastic Cuff	Yellow	110x106cm