



[Type 3B/4B Chemical Protective Gown]

Protection against numerous concentrated inorganic chemicals and biohazards. Protective seams, stitched and over-taped with barrier-tape, providing barrier performance equal to that of the fabric.

Fabric: Multi-layer laminated fabric

Seam Type: Stitched & taped seams

Style: One piece opening at back with ties, with tear tape closure on the neck, bottom sleeve with knitted cuffs (inner layer) and elastic cuffs (outer layer)

Standard:

EN 14605:2005+A1:2009 Type PB-3B/4B

Protective performance against spray and infective agents

EN 13034:2005 +A1:2009

Protective clothing against liquid chemicals

EN 14126:2003+AC:2004

Performance requirements and tests methods for protective clothing against infective agents

EN 1149-5:2018

Protective clothing – Electrostatic properties

EN 14325:2004

Protective clothing against chemicals - Test methods and performance classification of chemical protective clothing materials, seams, joins and assemblages

EN ISO 13688:2013 General requirements



REF	Color	Size
WLG1015	Yellow	S-3XL



EN 13034/05+A1/09
 Type [6]B
 EN 14605/05+A1/09
 Type [3]B / [4]B



EN 1149-5:2018



EN 14126:2003