



## [SMS Isolation Gown] WLG1002

SMS fabric is well-known for its balance between excellent water repellency and breathability. Modern workwear is one of the major playground for this fabric, especially used for cooling body from abundance of hot moisture.

Fabric: SMS 25-60gsm

Seam Type: Stitched/Ultrasonic welded

Style: One piece opening at back, ties at waist

With ties/velcro at back neck

Knitted/Elastic cuffs

Color: ● ○ ● ●

| REF        | Weight   | Style        | Color  | Size  |
|------------|----------|--------------|--------|-------|
| WLG1002    | 25-60gsm | Elastic Cuff | White  | S-3XL |
| WLG1002B   | 25-60gsm | Elastic Cuff | Blue   | S-3XL |
| WLG1002KC  | 25-60gsm | Knitted Cuff | White  | S-3XL |
| WLG1002KCB | 25-60gsm | Knitted Cuff | Blue   | S-3XL |
| WLG1002KCY | 25-60gsm | Knitted Cuff | Yellow | S-3XL |
| WLG1002Y   | 25-60gsm | Elastic Cuff | Yellow | S-3XL |

